




Gas Chromatography - Vacuum Ultraviolet Spectroscopy: A Versatile Tool for Analysis of Gasoline and Jet Fuels

Alex Hodgson, Jack Cochran, James Diekmann, Ryan Schonert

VUV Analytics, Inc.

PEFTEC 2019

Overview of Crude Oil Refining Products




- Liquefied petroleum gases (LPG): $C_1 - C_4$
- Gasoline: $C_5 - C_{12}$
- Middle distillates (e.g., jet, diesel): $C_{10} - C_{20}$
- Heavier oils and waxes: C_{20}^+

What is PIONA?

- **P**araffins (linear alkanes)




- **I**soparaffins (branched alkanes)




- **O**lefins (alkenes)



- **N**aphthenes (cycloalkanes)



- **A**romatics



Classical GC Detectors

Flame ionization detection (FID)

- Good LODs
- Wide dynamic range
- Cheap
- Robust
- Great for process applications

Mass spectrometry (MS)

- The “gold standard” (R&D)
- Compound identification
- Tunable physical selectivity (SIM, SRM, MRM)
- Very low LODs (MS/MS)


Shortcomings, and how VUV fits in



- FID requires baseline separation of peaks for accurate quantitation
- MS has limited flow rates due to vacuum pump considerations; can struggle with isomer differentiation
- ***VUV combines the speeds of FID with the qualitative spectral information of MS, and unique spectral fingerprints make deconvolution of coeluting isomers straightforward***

Vacuum Ultraviolet (VUV) Absorption Spectroscopy

- Absorption spectroscopy in a new spectral region (<200 nm)
 - Previously only measurable using synchrotrons
- A universal technique that provides unique spectral fingerprints
 - High energy, low wavelength exposure produces electronic transitions between $\sigma \rightarrow \sigma^*$, $n \rightarrow \sigma^*$, and $\pi \rightarrow \pi^*$ molecular orbitals
- ***Unambiguous compound identification, even for structural isomers***
- 1st order quantitation following Beer-Lambert Law ($A = \epsilon bc$)
- ***Chromatography runs can be deliberately compressed, leading to higher sample throughput***



	VGA-100	VGA-101
Max Temp	300°C	430°C
λ Range	125 - 240 nm	125 - 430 nm
Acq Rate	75 Hz	75 Hz



Traditional Gasoline Analysis

The current playing field of PIONA

DHA – Overview and Setup



- GC-FID (ASTM D6730), though sometimes simultaneously run on GC-MS for confirmation
- Speciation of up to 600 compounds in gasoline-range fuels (though not all compounds are identified by name)
- 100m non-polar column, 5% diphenyl “tuning” precolumn
- Cryogenic oven start (5°C)
- Run time: 174 minutes
 - Run time subsequently cut to 38 minutes

DHA – Drawbacks




- Relies solely on retention time
- Tuning precolumn necessary to baseline separate known coelutions
 - Benzene/1-Methylcyclopentene
 - m-/p-Xylene
- Any unexpected coelutions cannot be identified with the FID data, requires MS analysis
- Requires a trained user to go through data peak by peak

Reformulyzer – Overview and Setup


- Multi-dimensional GC-FID (ASTM D6839)
- Measures PONA compounds and oxygenates in gasoline and gasoline blend streams
 - Physically separates compounds using a series of columns and traps actuated by 6-port valves
 - Paraffins and isoparaffins not separated
- 39-minute run time




A "basic" schematic of the Reformulyzer...




A "basic" schematic of the Reformulyzer...




A "basic" schematic of the Reformulyzer...



A "basic" schematic of the Reformulyzer...



A "basic" schematic of the Reformulyzer...



Reformulyzer – Drawbacks

- Multiple columns and connections → greater chance for leaks or restrictions
- Access to certain areas only possible by removing all traps and fans, disconnecting columns
- Problems take longer to troubleshoot and repair
- Longer down time





PIONA by VUV: ASTM D8071

GC-VUV Analysis of Gasoline-Range Fuels



ASTM D8071-17

Standard Test Method for Determination of Hydrocarbon Group Types and Select Hydrocarbon and Oxygenate Compounds in Automotive Spark-Ignition Engine Fuel Using Gas Chromatography with Vacuum Ultraviolet Absorption Spectroscopy Detection (GC-VUV)

Active Standard ASTM D8071 | Developed by Subcommittee: [D02.04.01](#)

Paraffins, Isoparaffins, Olefins, Naphthenes, Aromatics by Carbon Number from C1 to C15

Methanol, Ethanol, Benzene, Toluene, Ethylbenzene, Xylenes, Naphthalene, Methylnaphthalenes, Isooctane

Recommended

[Standards Tracker](#)

[Standards Subscriptions](#)

Standards & Publications

All Standards & Publications

Standards Products

Symposia Papers & STPs

Manuals, Monographs, & Data Series

Journals

Reading Room

Authors

Book of Standards

Reading Room

Product Updates

Catalogs

Digital Library



Enterprise Solutions






PIONA Analysis of Gasoline with GC-VUV




- 30m x 0.25mm x 0.25 μ m 100% dimethylpolysiloxane column
 - About 33.5 min run versus two hours for DHA
 - No precolumn for benzene/1-methylcyclopentene separation, e.g.
- Collect VUV absorbance data (125-240 nm)
 - Spectra can be unique and class-indicating
- Fully automated data processing with VUV Analyze PIONA +
 - Spectral library, retention index, relative response factors
 - "Separation" of coelutions with Time Interval Deconvolution
 - Mass % and Volume % reports, carbon number breakdown

Precolumn? We don't need no stinkin' precolumn...








Olefin/Naphthene Differentiation using VUV Spectroscopy



Issues with Olefin/Naphthene Differentiation using MS




ASTM D8071 Run Conditions


- Agilent 6890 GC

- 1 μL injection
- Inlet: 250°C, split 300:1
- Column: Restek Rxi-1ms (30 m x 0.25 mm, 0.25- μm)
- Flow: 1 mL/min helium, constant flow
- Oven: 35°C, hold 10 min; 7°C/min to 200°C (run time – 33.6 min)

- VUV Analytics VGA-100

- Makeup gas: 0.40 psi N_2
- Flow cell and transfer line: 275°C
- Acquisition rate: 4.5 spectra/sec
- Acquisition range: 125-240 nm





Click to add or remove plots

**Changes to the legend are not persistent

125 - 240 (8071)



R & D Results

Report File

View Report

Export to Excel

Input Files

Reverse Search

Create RI File

D8071

Analytes to Include

Advanced

Spectra

Fit Info

Run File (*.db) or Directory C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\Demo VUV-CS.db



Refresh

Retention Index Markers File C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\Demo VUV RI.txt



Refresh

Reference Library File C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\WUV PIONA ASTM D8071 rev 1.1.0.db



Refresh

Report Method C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\D8071-18 Report Method.VUV_RM



Refresh

Edit / Create ...

Analysis Parameters

Initial Background Time (min)

Begin 1.600

End 1.800

Analysis Time (min)

Begin 1.800

End 32.000

Time Step 0.020

Methods

PIONA

Tiered Search Limit

3 Analytes

Chromatogram Filter

8071 Filters

Chi² Min

1.0000E-9

Chi² Max

1.0000E-1

Chi² delta (%)

60.00

R² Limit

0.4000

Abs Threshold

0.00050

BG Threshold

0.00025

BG Scalar

1.5

RI window +/-

25


Analyze

Load Parameters

Time Step Status

Save Parameters

Stop Analyzing



Click to add or remove plots

**Changes to the legend are not persistent

125 - 240 (8071)

R & D Results Report File View Report Export to Excel

File from gasoline analysis with GC-VUV

Input Files

Run File (*.db) or Directory

Retention Index Markers File C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\Demo VUV RI.txt


Reference Library File C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\WUV PIONA ASTM D8071 rev 1.1.0.db

Report Method C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\D8071-18 Report Method.VUV_RM

Analysis Parameters

Initial Background Time (min)	Analysis Time (min)	Methods	Tiered Search Limit	Chi ² Min	Chi ² delta (%)	Abs Threshold	BG Scalar
Begin 1.600	Begin 1.800	PIONA	3 Analytes	1.0000E-9	60.00	0.00050	1.5
End 1.800	End 32.000		Chromatogram Filter	Chi ² Max	R ² Limit	BG Threshold	RI window +/-
Time Step 0.020			8071 Filters	1.0000E-1	0.4000	0.00025	25

Analyze Load Parameters Time Step Status Save Parameters Stop Analyzing



Click to add or remove plots

**Changes to the legend are not persistent

125 - 240 (8071)

R & D Results Report File View Report Export to Excel

Input Files Reverse Search Create RI File D8071 Analytes to Include Advanced Spectra Fit Info

Run File (*.db) or Directory Refresh

Retention Index Markers File **Retention Index markers file** Refresh

Reference Library File C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\WUV PIONA ASTM D8071 rev 1.1.0.db Refresh

Report Method C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\D8071-18 Report Method.VUV_RM Refresh Edit / Create ...

Analysis Parameters

Initial Background Time (min)

Begin 1.600 End 1.800 Time Step 0.020

Analysis Time (min)

Begin 1.800 End 32.000

Methods PIONA

Tiered Search Limit 3 Analytes

Chromatogram Filter 8071 Filters

Chi² Min 1.0000E-9

Chi² Max 1.0000E-1

Chi² delta (%) 60.00

R² Limit 0.4000


Abs Threshold 0.00050

BG Threshold 0.00025

BG Scalar 1.5

RI window +/- 25

Analyze Load Parameters Time Step Status Save Parameters Stop Analyzing



Click to add or remove plots

**Changes to the legend are not persistent

125 - 240 (8071)

R & D Results Report File View Report Export to Excel

Input Files Reverse Search Create RI File D8071 Analytes to Include Advanced Spectra Fit Info

Run File (*.db) or Directory C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\Demo VUV-CS.db Refresh

Retention Index Markers File Refresh


Reference Library File **PIONA compound absorbance spectra library**

Report Method C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\18 Report Method.VUV_RM Refresh Edit / Create ...

Analysis Parameters

Initial Background Time (min)	Analysis Time (min)	Methods	Tiered Search Limit	Chi ² Min	Chi ² delta (%)	Abs Threshold	BG Scalar
Begin 1.600	Begin 1.800	PIONA	3 Analytes	1.0000E-9	60.00	0.00050	1.5
End 1.800	End 32.000		Chromatogram Filter	Chi ² Max	R ² Limit	BG Threshold	RI window +/-
Time Step 0.020			8071 Filters	1.0000E-1	0.4000	0.00025	25

Analyze Load Parameters Time Step Status Save Parameters Stop Analyzing



Click to add or remove plots

**Changes to the legend are not persistent

125 - 240 (8071)

R & D Results Report File View Report Export to Excel

Input Files Reverse Search Create RI File D8071 Analytes to Include Advanced Spectra Fit Info

Run File (*.db) or Directory C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\Demo VUV-CS.db Refresh

Retention Index Markers File C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\Demo VUV RI.txt Refresh


Reference Library File C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\WUV PIONA ASTM D8071 rev 1.1.0.db Refresh

Report Method **PIONA report method** V_RM Refresh Edit / Create ...

Analysis Parameters

Initial Background Time (min)	Analysis Time (min)	Methods	Tiered Search Limit	Chi ² Min	Chi ² delta (%)	Abs Threshold	BG Scalar
Begin 1.600	Begin 1.800	PIONA	3 Analytes	1.0000E-9	60.00	0.00050	1.5
End 1.800	End 32.000		Chromatogram Filter	Chi ² Max	R ² Limit	BG Threshold	RI window +/-
Time Step 0.020			8071 Filters	1.0000E-1	0.4000	0.00025	25

Analyze Load Parameters Time Step Status Save Parameters Stop Analyzing



Click to add or remove plots
**Changes to the legend are not persistent
 125 - 240 (8071)

Single File or Batch Processing

R & D Results Report File View Report Export to Excel

Input Files	Reverse Search	Create RI File	D8071	Analytes to Include	Advanced	Spectra	Fit Info
<input type="text" value="Analytics\WUV Analyze 130 1953\Demos\D8071\Demo VUV-CS.db"/>						<input type="button" value="Refresh"/>	
<input type="text" value="Analytics\WUV Analyze 130 1953\Demos\D8071\Demo VUV RI.txt"/>						<input type="button" value="Refresh"/>	
<input type="text" value="Analytics\WUV Analyze 130 1953\Demos\WUV PIONA ASTM D8071 rev 1.1.0.db"/>						<input type="button" value="Refresh"/>	
<input type="text" value="Report Method C:\Program Files\WUV Analytics\WUV Analyze 130 1953\Demos\D8071\18 Report Method.VUV_RM"/>						<input type="button" value="Refresh"/> <input type="button" value="Edit / Create ..."/>	

Analysis Parameters


Initial Background (min)	Analysis Time (min)	Methods	Tiered Search Limit	Chi ² Min	Chi ² delta (%)	Abs Threshold	BG Scalar
Begin 1.800	Begin 1.800	PIONA	3 Analytes	1.0000E-9	60.00	0.00050	1.5
End 1.800	End 32.000		Chromatogram Filter	Chi ² Max	R ² Limit	BG Threshold	RI window +/-
Time Step 0.020			8071 Filters	1.0000E-1	0.4000	0.00025	25

Time Interval Deconvolution (TID)



- Alternative quantitation method using VUV Analyze™ software
- Chromatogram is divided into equal, small time intervals (typically <0.05 min)
- For each time interval, compare measured spectrum against reference spectra in designated library, best analyte(s) fit determined
- Can quickly be performed to measure total response per analyte for a chromatogram; this can be converted into a relative mass percentage and relative volume percentage

VUV Spectral Deconvolution of Hydrocarbons in Gasoline




Chi² 2: 1.735E-4


R²: 0.9992

Analyte Table

RRF Display	Analyte Name	Analyte Category	Fit Values	Spectra Ave.	Response	Area	Ret. Index
Olefin	trans-2-Pentene	Olefin	0.85770	0.35124	0.30126	0.001115	506.73599
Paraffin	Pentane	Linear alkane	0.96550	0.18176	0.17549	0.0006493	500.00000
	2-Methyl-1,3-butadiene	Olefin	0.36569	0.45053	0.16476	0.0006096	503.70239

Spectral Deconvolution of Coeluting Hydrocarbons





Click to add or remove plots

**Changes to the legend are not persistent

- 125 - 240 (8071)
- Paraffin
- Isoparaffin
- Olefin
- Naphthene
- Aromatic
- Methanol
- Ethanol
- Methyl tert-butyl ether
- iso-octane
- Naphthalene
- Methylnaphthalenes
- Benzene
- Toluene
- Ethylbenzene

Mass % Report for PIONA + Select Individual Hydrocarbons and Ethanol


Response	Mass %					Volume %	Amount				Config	
	P	I	O	N	A		Category	Mass %	RT (min.)	RRF		Density
C1							Paraffin	8.8192				
C2							Isoparaffin	42.5362				
C3							Olefin	9.1821				
C4	1.2604	0.1114	0.0581				Naphthene	8.1559				
C5	3.1375	6.7534	3.7404	0.1456			Aromatic	20.4855				
C6	1.9858	7.3473	2.8238	1.6022	0.8250		Methanol	0.0000				
C7	1.1744	5.6581	1.0126	2.2271	3.9218		Ethanol	10.8211	2.370	1.029	0.789	
C8	0.5486	18.1335	1.0130	1.6267	5.7333		Methyl tert-butyl ether	0.0000				
C9	0.3854	2.8188	0.1990	1.3372	5.8269		iso-octane	7.5817	4.930	0.674	0.692	
C10	0.1992	0.8996	0.2441	0.7377	2.6590		Naphthalene	0.1886	22.429	0.207	1.025	
C11	0.0894	0.6189	0.0679	0.4207	1.2177		Methylnaphthalenes	0.2171	24.810	0.250	1, 1.010	
C12	0.0385	0.1390	0.0232	0.0588	0.2667		Benzene	0.8250	4.029	0.258	0.879	
C13		0.0514			0.0350		Toluene	3.9218	7.209	0.267	0.867	
C14		0.0047					Ethylbenzene	0.9860	12.730	0.284	0.867	
C15							Xylenes	4.7461	13.190	0.284	0.870	
C16												
C17												
C18												
C19												
Total	8.8192	42.5362	9.1821	8.1559	20.4855							

Area Rejected % (Chi² Limit)

Negative Fit %

No. of Background Time steps

Total No. of Time steps



Click to add or remove plots

**Changes to the legend are not persistent

- 125 - 240 (8071)
- Paraffin
- Isoparaffin
- Olefin
- Naphthene
- Aromatic
- Methanol

Response			Mass %			Volume %		Amount			Config	
P	I	O	N	A			Category	Mass %	RT (min.)	RRF	Density	
C1							Paraffin	8.8192				
C2							Isoparaffin	42.5362				
C3							Olefin	9.1821				
C4	1.2604	0.1114	0.0581				Naphthene	8.1559				
C5	3.1375	6.7534	3.7404	0.1456			Aromatic	20.4855				
C6	1.9858	7.3473	2.8238	1.6022	0.8250		Methanol	0.0000				
C7	1.1744	5.6581	1.0126	2.2271	3.9218		Ethanol	10.8211	2.370	1.029	0.789	
C8	0.5486	18.1335	1.0130	1.6267	5.7333		Methyl tert-butyl ether	0.0000				
C9	0.3854	2.8188	0.1990	1.3372	5.8269		iso-octane	7.5817	4.930	0.674	0.692	
C10	0.1992	0.8996	0.2441	0.7377	2.6590		Naphthalene	0.1886	22.429	0.207	1.025	
C11	0.0894	0.6189	0.0679	0.4207	1.2177		Methylnaphthalenes	0.2171	24.810	0.250	1, 1.010	
C12	0.0385	0.1390	0.0232	0.0588	0.2667		Benzene	0.8250	4.029	0.258	0.879	
C13		0.0514			0.0350		Toluene	3.9218	7.209	0.267	0.867	
C14		0.0047					Ethylbenzene	0.9860	12.730	0.284	0.867	
C15							Xylenes	4.7461	13.190	0.284	0.870	
C16												
C17												
C18												
C19												
Total	8.8192	42.5362	9.1821	8.1559	20.4855							

List Poor Fits

C17					
C18					
C19					
Total	8.8192	42.5362	9.1821	8.1559	20.4855




Conjugated Diolefin Analysis by GC-VUV

Putting the “+” in PIONA+™

The Trouble with Conjugated Diolefins (CDOs)

- Easily polymerize in high enough concentrations
 - Used in synthesis of rubber and other polymers
- Polymers can foul hydrocarbon streams of petrochemical refineries
- CDO levels must be kept below a threshold to prevent polymerization

1,2-Polymerization of 1,3- Butadiene




1,4-Polymerization of 1,3- Butadiene



UOP-326: Maleic Anhydride Method


- Common method for determining diene content in lighter hydrocarbon mixtures
- Maleic anhydride consumed via Diels-Alder reaction with dienes; remaining maleic anhydride converted to maleic acid, which is measured by colorimetric titration



Drawbacks of Using Maleic Anhydride

- >3-hour analysis time (both manual and automated)
- Nucleophiles (e.g. alcohols, thiols) will also react with maleic anhydride; high diene values reported for samples containing MeOH and EtOH
- Some sterically-hindered dienes (e.g. 2,5-dimethyl-2,4-hexadiene) will not react
- Only semiquantitative, as reaction is not selective
- No qualitative information (i.e. identification of diolefin species)

VUV Absorbance Spectra for Conjugated Diolefins



Conjugated Diolefin Run Conditions

- Agilent 6890 GC

- 1 μL injection
- Inlet: 250°C, split 300:1
- Column: Restek Rxi-1ms (30 m x 0.25 mm, 0.25- μm)
- Flow: 1 mL/min helium, constant flow
- Oven: 35°C, hold 10 min; 7°C/min to 200°C (run time – 33.6 min)

- VUV Analytics VGA-100

- Makeup gas: 0.40 psi N_2
- Flow cell and transfer line: 275°C
- Acquisition rate: 4.5 spectra/sec
- Acquisition range: 125-240 nm




Conjugated Diolefin Run Conditions

- Agilent 6890 GC
 - 1 μ L injection
 - Inlet: 250°C, split 300:1
 - Column: Restek Rxi-1ms (5) m \times 0.25 mm, 0.25- μ m
 - Flow: 1 mL/min helium constant flow
 - Oven: 35°C hold 10 min; 7°C/min to 210°C (run time – 33.6 min)
- VUV Analytical VGA-100
 - Makeup gas: 0.40 psi N₂
 - Flow cell and transfer line: 275°C
 - Acquisition rate: 4.5 spectra/sec
 - Acquisition range: 125-240 nm



**SAME METHOD AS
ASTM D8071!**



Gasoline Spiked with 0-1% Conjugated Diolefins



Spiked Gasoline Overlays: 0-1% CDOs




Quantitative Analysis of CDO-Spiked Gasoline

2-Methyl-1,3-butadiene		1,3-Cyclohexadiene		2,4-Dimethyl-1,3-pentadiene		2,5-Dimethyl-2,4-hexadiene		1,3-Cyclooctadiene	
Expected Conc	Measured Conc	Expected Conc	Measured Conc	Expected Conc	Measured Conc	Expected Conc	Measured Conc	Expected Conc	Measured Conc
0	0.015	0	-	0	-	0	-	0	-
0.010	0.022	0.010	-	0.010	-	0.010	-	0.010	-
0.020	0.029	0.020	0.019	0.020	-	0.020	-	0.020	-
0.050	0.047	0.051	0.044	0.050	0.042	0.051	0.048	0.050	0.028
0.10	0.11	0.10	0.12	0.10	0.10	0.10	0.11	0.10	0.099
0.20	0.21	0.20	0.23	0.20	0.22	0.20	0.22	0.20	0.24
0.50	0.48	0.51	0.56	0.50	0.55	0.51	0.52	0.50	0.58
1.0	1.1	1.0	1.2	1.0	1.2	1.0	1.1	1.0	1.3

*All concentrations in mass %

Linearity for Conjugated Diolefins in Gasoline





VUV Verified™ Hydrocarbon Analysis (VHA™)

Detailed Hydrocarbon Analysis with GC-VUV

The Starting Point for VUV VHA – ASTM D6730



ASTM INTERNATIONAL
Helping our world work better

All



Search topic, title, author, A53



SIGN IN



PRODUCTS & SERVICES | GET INVOLVED | ABOUT | HOME | NEWS

Languages | Contact | Cart

Products and Services / Standards & Publications / Standards Products

Standards & Publications

All Standards and
Publications

Standards Products

Symposia Papers & STPs

Manuals, Monographs, &
Data Series

Journals

ASTM D6730 - 01(2016) ⓘ

Standard Test Method for Determination of Individual
Components in Spark Ignition Engine Fuels by 100–Metre
Capillary (with Precolumn) High-Resolution Gas Chromatography

Active Standard ASTM D6730 | Developed by Subcommittee: [D02.04.01](#)

Recommended

Standards Tracker

Standards Subscriptions

Customers frequently
also buy:


D5453

Strategy for Prototype DHA (VHA) with GC-VUV

- Use ASTM D6730 GC conditions
 - 100m x 0.25mm x 0.50 μ m “methyl silicone” GC column
 - Exception; **no precolumn** for critical separations
 - Cryogenic oven start (5°C)
- Build VUV VHA gasoline compound retention time table using D6730 DHA Component Table
 - Starting with higher mass % compounds in gasoline
- Use VUV Analyze software for automated data processing
 - “Analytes to Include”
 - Time Interval Deconvolution

GC Conditions – VUV VHA

- Agilent 6890 GC
- Precision 4mm Split liner with Wool (Restek)
- 250°C, 1 μ L
- Split ratio 300:1
- 100m x 0.25mm x 0.50 μ m Rtx-DHA GC Column (Restek)
- Helium 40 psi (constant pressure)
- 5°C (10 min), 5°C/min to 50°C (50 min), 1.5°C/min to 200°C (5 min)
- **Run time: 174 min**



VUV Spectrometer Conditions – VUV VHA


- VUV Analytics VGA-100
- Transfer line temperature: 275°C
- Flow cell temperature: 275°C
- Makeup gas pressure: 0.25 psi
- Acquisition range: 125 to 240 nm
- Acquisition rate: 4 spectra/sec



Status Reference Library: C:\Users\Jack Cochran\Documents\VUV Data\D6730 20180224\VUV Analyze Files D6730\VUV Iib D6730 20180224.db
Last Mod Date: 2018/02/23 14:16:24

Current Run File C:\Users\Jack Cochran\Documents\VUV Data\D6730 20180224\PONA VI and VIIData Files\VUV1 D6730 0223 054 08 25 22.db

Production Analyze Results



125 - 240 Spectra Ave.
 125 - 240 Avg.
 Aromatics
 Saturates
 Di-Olefins

Fit Info 0.000
 Background Begin 3.000
 Background End 4.000
 Analyze Begin 5.000
 Analyze End 170.0

Input Files Reverse Search Create RI File **Analytes to Include** Advanced Peaks Spectra Fit Info

Name	Include in the...	Start time (min)	End time (min)	Use during
Butane	List of possible analytes	9.22	9.62	Only during time range
iso-pentane	List of possible analytes	12.83	13.23	Only during time range
1-Pentene	List of possible analytes	13.81	14.21	Only during time range
2-Methyl-1-butene	List of possible analytes	14.28	14.68	Only during time range
Pentane	List of possible analytes	14.57	14.97	Only during time range
trans-2-Pentene	List of possible analytes	15.14	15.54	Only during time range

Analysis Parameters


Initial Background Time (min) Begin End
 Analysis Time (min) Begin End Time Step
 Methods PIONA Oxygenates
 DHA
 Tiered Search Limit 3 Analytes Chromatogram Filter 140 - 160
 Chi² Min 1.0000E-9 Chi² Max 1.0000E-1
 Chi² delta (%) 40 R² Limit 0.8000
 Abs Threshold 0.0010 BG Threshold 0.0003
 Use Peak Detection
 Analyze Spectra Within Peaks
 Use Initial Background Time
 RI window +/-

Analyze Load Parameters Save Parameters Stop Analyzing

Status Reference Library: C:\Users\Jack Cochran\Documents\VUV Data\D6730 20180224\VUV Analyze Files D6730\VUV Iib D6730 20180224.db
Last Mod Date: 2018/02/23 14:16:24

Current Run File C:\Users\Jack Cochran\Documents\VUV Data\D6730 20180224\PONA VI and VIIData Files\VUV1 D6730 0223 054 08 25 22.db

Production Analyze Results



125 - 240 Spectra Ave.
 125 - 240 Avg.
 Aromatics
 Saturates
 Di-Olefins

Fit Info 0.000
 Background Begin 3.000
 Background End 4.000
 Analyze Begin 5.000
 Analyze End 170.0

Input Files Reverse Search Create RI File **Analytes to Include** Advanced Peaks Spectra Fit Info




Name	Include in the...	Start time (min)	End time (min)	Use during
1,2-Dimethyl-3-ethylbenzene	List of possible analytes	121.22	121.62	Only during time range
1,2,4,5-Tetramethylbenzene	List of possible analytes	122.91	123.31	Only during time range
1,2,3,5-Tetramethylbenzene	List of possible analytes	123.36	123.76	Only during time range
Naphthalene	List of possible analytes	130.49	130.89	Only during time range
2-Methylnaphthalene	List of possible analytes	143.47	143.87	Only during time range
1-Methylnaphthalene	List of possible analytes	145.15	145.55	Only during time range





Analysis Parameters

Initial Background Time (min) Begin End
 Analysis Time (min) Begin End Time Step
 Methods PIONA Oxygenates
 DHA
 Tiered Search Limit 3 Analytes Chromatogram Filter 140 - 160
 Chi² Min 1.0000E-9 Chi² Max 1.0000E-1
 Chi² delta (%) 40 R² Limit 0.8000
 Abs Threshold 0.0010 BG Threshold 0.0003
 Use Peak Detection
 Analyze Spectra Within Peaks
 Use Initial Background Time
 RI window +/-

Analyze Load Parameters Save Parameters Stop Analyzing

VUV VHA
No precolumn











VUV VHA - By the Numbers



- Currently report mass and volume % data for 151 compounds
- Individual compounds add up to about 90% of gasoline composition
 - Other compounds are classed through absorbance spectra shape and RI
 - Class carbon number report from C₃ to C₁₅
- VHA (PIONA) absorbance spectra library contains ~ 750 spectra
 - Additional spectra always being added

Class	#
Paraffins	13
Isoparaffins	27
Olefins	45
Naphthenes	24
Aromatics	28
Di-Aromatics	9
Alcohols	2
Ethers	3









Mass %

Category	MEG D8071	MEG VHA		B02 D8071	B02 VHA		M872 D8071	M872 VHA		M873 D8071	M873 VHA
Paraffin	13.7	15.2		17.6	17.5		7.93	8.18		12.8	12.9
Isoparaffin	27.6	28.5		18.9	18.7		68.6	70.0		25.1	24.5
Olefin	13.2	12.7		13.4	14.5		1.98	0.929		0.276	0.003
Naphthene	8.84	6.37		2.26	1.84		2.12	1.06		5.65	5.54
Aromatic	25.6	24.8		36.6	36.1		19.3	19.9		55.4	56.3
Di-Aromatic	0.970	0.859		0.063			0.152			0.788	0.760
Ethanol	10.1	11.5		11.1	11.4						
Iso-octane	1.91	2.08		4.87	5.07		16.6	17.7			
Benzene	0.426	0.479					0.222	0.234		0.446	0.451
Toluene	3.51	3.99		15.4	16.0		6.91	7.70		15.1	15.5
Ethylbenzene	0.875	0.847		1.37	1.15		1.37	1.31		3.25	3.03
Xylenes	5.55	5.83		6.50	6.81		6.88	7.40		16.1	17.2
Naphthalene	0.266	0.279		0.063			0.053			0.401	0.437
Methylnaphthalenes	0.590	0.580					0.099				

Mass %

Category	MEG D8071	MEG VHA		B02 D8071	B02 VHA		M872 D8071	M872 VHA		M873 D8071	M873 VHA
Paraffin	13.7	15.2		17.6	17.5		7.93	8.18		12.8	12.9
Isoparaffin	27.6	28.5		18.9	18.7		68.6	70.0		25.1	24.5
Olefin	13.2	12.7		13.4	14.5		1.98	0.929		0.276	0.003
Naphthene	8.84	6.37		2.26	1.84		2.12	1.06		5.65	5.54
Aromatic	25.6	24.8		36.6	36.1		19.3	19.9		55.4	56.3
Di-Aromatic	0.970	0.859		0.063			0.152			0.788	0.760
Ethanol	10.1	11.5		11.1	11.4						
Iso-octane	1.91	2.08		4.87	5.07		16.6	17.7			
Benzene	0.426	0.479					0.222	0.234		0.446	0.451
Toluene	3.51	3.99		15.4	16.0		6.91	7.70		15.1	15.5
Ethylbenzene	0.875	0.847		1.37	1.15		1.37	1.31		3.25	3.03
Xylenes	5.55	5.83		6.50	6.81		6.88	7.40		16.1	17.2
Naphthalene	0.266	0.279		0.063			0.053			0.401	0.437
Methylnaphthalenes	0.590	0.580					0.099				

Mass %

Category	MEG D8071	MEG VHA		B02 D8071	B02 VHA		M872 D8071	M872 VHA		M873 D8071	M873 VHA
Paraffin	13.7	15.2		17.6	17.5		7.93	8.18		12.8	12.9
Isoparaffin	27.6	28.5		18.9	18.7		68.6	70.0		25.1	24.5
Olefin	13.2	12.7		13.4	14.5		1.98	0.929		0.276	0.003
Naphthene	8.84	6.37		2.26	1.84		2.12	1.06		5.65	5.54
Aromatic	25.6	24.8		36.6	36.1		19.3	19.9		55.4	56.3
Di-Aromatic	0.970	0.859		0.063			0.152			0.788	0.760
Ethanol	10.1	11.5		11.1	11.4						
Iso-octane	1.91	2.08		4.87	5.07		16.6	17.7			
Benzene	0.426	0.479					0.222	0.234		0.446	0.451
Toluene	3.51	3.99		15.4	16.0		6.91	7.70		15.1	15.5
Ethylbenzene	0.875	0.847		1.37	1.15		1.37	1.31		3.25	3.03
Xylenes	5.55	5.83		6.50	6.81		6.88	7.40		16.1	17.2
Naphthalene	0.266	0.279		0.063			0.053			0.401	0.437
Methylnaphthalenes	0.590	0.580					0.099				

Mass %

Category	MEG D8071	MEG VHA		B02 D8071	B02 VHA		M872 D8071	M872 VHA		M873 D8071	M873 VHA
Paraffin	13.7	15.2		17.6	17.5		7.93	8.18		12.8	12.9
Isoparaffin	27.6	28.5		18.9	18.7		68.6	70.0		25.1	24.5
Olefin	13.2	12.7		13.4	14.5		1.98	0.929		0.276	0.003
Naphthene	8.84	6.37		2.26	1.84		2.12	1.06		5.65	5.54
Aromatic	25.6	24.8		36.6	36.1		19.3	19.9		55.4	56.3
Di-Aromatic	0.970	0.859		0.063			0.152			0.788	0.760
Ethanol	10.1	11.5		11.1	11.4						
Iso-octane	1.91	2.08		4.87	5.07		16.6	17.7			
Benzene	0.426	0.479					0.222	0.234		0.446	0.451
Toluene	3.51	3.99		15.4	16.0		6.91	7.70		15.1	15.5
Ethylbenzene	0.875	0.847		1.37	1.15		1.37	1.31		3.25	3.03
Xylenes	5.55	5.83		6.50	6.81		6.88	7.40		16.1	17.2
Naphthalene	0.266	0.279		0.063			0.053			0.401	0.437
Methylnaphthalenes	0.590	0.580					0.099				

Mass %

Category	MEG D8071	MEG VHA		B02 D8071	B02 VHA		M872 D8071	M872 VHA		M873 D8071	M873 VHA
Paraffin	13.7	15.2		17.6	17.5		7.93	8.18		12.8	12.9
Isoparaffin	27.6	28.5		18.9	18.7		68.6	70.0		25.1	24.5
Olefin	13.2	12.7		13.4	14.5		1.98	0.929		0.276	0.003
Naphthene	8.84	6.37		2.26	1.84		2.12	1.06		5.65	5.54
Aromatic	25.6	24.8		36.6	36.1		19.3	19.9		55.4	56.3
Di-Aromatic	0.970	0.859		0.063			0.152			0.788	0.760
Ethanol	10.1	11.5		11.1	11.4						
Iso-octane	1.91	2.08		4.87	5.07		16.6	17.7			
Benzene	0.426	0.479					0.222	0.234		0.446	0.451
Toluene	3.51	3.99		15.4	16.0		6.91	7.70		15.1	15.5
Ethylbenzene	0.875	0.847		1.37	1.15		1.37	1.31		3.25	3.03
Xylenes	5.55	5.83		6.50	6.81		6.88	7.40		16.1	17.2
Naphthalene	0.266	0.279		0.063			0.053			0.401	0.437
Methylnaphthalenes	0.590	0.580					0.099				

	P	I	O	N	A
C0					
C1					
C2					
C3					
C4	5.1914	0.2350	0.6322		
C5	3.5075	7.5802	6.8148	0.1822	
C6	2.1806	6.6897	3.7730	1.8755	0.4786
C7	1.5994	4.7750	0.9201	2.8040	3.9930
C8	1.0650	6.4443	0.5264	1.3166	6.6773
C9	0.7974	1.5549	0.0032	0.0670	9.1510
C10	0.4157	0.9806	0.0298	0.0140	3.5778
C11	0.2058	0.2394	0.0094	0.1076	0.8498
C12	0.1420	0.0058			0.0759
C13	0.0784	0.0115			0.0460
C14					
C15					
Total	15.1833	28.5164	12.7089	6.3669	24.8493

Category	Mass %
Paraffin	15.1833
Isoparaffin	28.5163
Olefin	12.7089
Naphthene	6.3669
Aromatic	24.8493
Di-Aromatic	0.8587
iso-butane	0.2350
Methanol	0.0000
Isobutylene	0.0268
1-Butene	0.0424
1,3-Butadiene	0.0000
Butane	5.1914
Vinyl acetylene	0.0000
trans-2-Butene	0.2722
cis-2-Butene	0.2909
Ethyl Alcohol	11.5165
3-Methyl-1-butene	0.1967
iso-pentane	7.5605
1,4-Pentadiene	0.0000
1-Pentene	0.5959
2-Methyl-1-butene	1.1212

Gasoline analyzed on 100m x 0.25mm x 0.50µm Rtx-DHA using D6730 GC conditions and VUV VHA data processing...

Category	Mass %
2,4-Dimethylpentane	0.4238
1-Methylcyclopentene	0.3283
Benzene	0.4786
Cyclohexane	0.3414
4-Methyl-1-hexene	0.0000
cis/trans-4-Methyl-2-hexene	0.1032
2-Methylhexane	1.9064
2,3-Dimethylpentane	0.4675
1,1-Dimethylcyclopentane	0.0329
tert-Amyl methyl ether	0.0000
Cyclohexene	0.0000
3-Methylhexane	1.7435
cis-1,3-Dimethylcyclopentane	0.5218
trans-1,3-Dimethylcyclopentane	0.4231
3-Ethylpentane	0.1286
trans-1,2-Dimethylcyclopentane	0.2884
iso-octane	2.0745
1-Heptene	0.0000
trans-3-Heptene	0.0247
Heptane	1.5994
trans-2-Heptene	0.1160

Category	Mass %
2,5-Dimethylhexane	0.3725
2,4-Dimethylhexane	0.4312
ctc-1,2,4-Trimethylcyclopentane	0.1041
ctc-1,2,3-Trimethylcyclopentane	0.0000
2,3,3-Trimethylpentane	0.2016
Toluene	3.9930
2,3-Dimethylhexane	0.2884
1,1,2-Trimethylcyclopentane	0.0000
2-Methylheptane	0.0000
4-Methylheptane	0.2470
3-Methylheptane	0.9112
trans-1,4-Dimethylcyclohexane	0.0991
2,2,5-Trimethylhexane	0.1434
3c-Ethylmethylcyclopentane	0.0744
3t-Ethylmethylcyclopentane	0.0900
2t-Ethylmethylcyclopentane	0.0000
trans-1,2-Dimethylcyclohexane	0.0171
Octane	1.0650
cis-1,2-Dimethylcyclohexane	0.0000
1,1,4-Trimethylcyclohexane	0.0000
Ethylcyclohexane	0.1275

Category	Mass %
2,5-Dimethylheptane	0.0473
Ethylbenzene	0.8470
ctt-1,2,4-Trimethylcyclohexane	0.0000
m-Xylene	2.9707
p-Xylene	1.1242
2-Methyloctane	0.4273
Styrene	0.0000
3-Methyloctane	0.3826
o-Xylene	1.7354
Nonane	0.7974
Isopropylbenzene	0.1543
n-Propylbenzene	0.4698
1-Methyl-3-ethylbenzene	1.9510
1-Methyl-4-ethylbenzene	0.7819
1,3,5-Trimethylbenzene	1.0287
1-Methyl-2-ethylbenzene	0.6209
2-Methylnonane	0.2587
3-Methylnonane	0.2387
1,2,4-Trimethylbenzene	3.4281
Decane	0.4157
1,2,3-Trimethylbenzene	0.7162

Category	Mass %
Dicyclopentadiene	0.0000
Indene	0.0000
1-Methyl-3-n-propylbenzene	0.3674
1-Methyl-4-n-propylbenzene	0.1724
n-Butylbenzene	0.0000
1,3-Dimethyl-5-ethylbenzene	0.4231
trans-Decalin	0.0000
1-Methyl-2-n-propylbenzene	0.0944
1,4-Dimethyl-2-ethylbenzene	0.2814
1,3-Dimethyl-4-ethylbenzene	0.2362
1,2-Dimethyl-4-ethylbenzene	0.4279
1,2-Dimethyl-3-ethylbenzene	0.0947
Undecane	0.2058
1,2,4,5-Tetramethylbenzene	0.2434
1,2,3,5-Tetramethylbenzene	0.3710
Naphthalene	0.2792
Dodecane	0.1420
Pentamethylbenzene	0.0149
2-Methylnaphthalene	0.4107
1-Methylnaphthalene	0.1688
Tridecane	0.0784

Can We Improve VUV VHA GC Run Time?

- Shorter GC column
- Thinner column film thickness
- Method translation for a new GC oven program
- Elution order of compounds should be the same
- GREATLY reduced run time
- What about the separations?

EZGC™ Method Translator

Carrier Gas	Original	Translation
	Helium	Helium

Column

Length	100.00	60.00 m
Inner Diameter	0.25	0.25 mm
Film Thickness	0.50	0.25 µm
Phase Ratio	125	250

Control Parameters

Outlet Flow	2.07	→ 2.07 mL/min
Average Velocity	25.50	31.27 cm/sec
Holdup Time	6.54	3.20 min
Inlet Pressure (gauge)	→ 40.00	28.52 psi
Outlet Pressure (abs)	14.70	14.70 psi

Atm Vacuum Atm Vacuum

Oven Program

Isothermal Ramps

	Ramp (°C/min)	Temp (°C)	Hold (min)	Ramp (°C/min)	Temp (°C)	Hold (min)
Number of Ramps		5	10		5	4.05
2 (1-4)	5	50	50	18.8	50	13.15
	1.5	200	5	5.4	200	1.4



Control Method



Constant Flow



Results



Solve for Efficiency Speed Translate Custom

Run Time	174.00	48.77 min
Speed		3.57 x









Mass % Comparison for Gasoline Analyzed with 3 GC-VUV Methods

- ASTM D8071 PIONA
 - 30m x 0.25mm x 0.25 μ m, no cryo
- ASTM D6730 VUV VHA
 - 100m x 0.25mm x 0.50 μ m, cryo
- "Fast" VUV VHA
 - 60m x 0.25mm x 0.25 μ m, cryo

Category	30m D8071	100m VHA	60m VHA
Paraffin	13.7	15.2	13.8
Isoparaffin	27.6	28.5	27.7
Olefin	13.2	12.7	12.3
Naphthene	8.84	6.37	7.75
Aromatic	25.6	24.8	26.9
Di-Aromatic	0.970	0.859	0.994
Ethanol	10.1	11.5	10.6
Iso-octane	1.91	2.08	2.00
Benzene	0.426	0.479	0.441
Toluene	3.51	3.99	3.82
Ethylbenzene	0.875	0.847	0.913
Xylenes	5.55	5.83	5.96
Naphthalene	0.266	0.279	0.345
Methylnaphthalenes	0.590	0.580	0.615



Jet Fuel Analysis by GC-VUV

Flying at the Speed of (Ultraviolet) Light

Jet Fuel Analysis – An Overview

- AKA Aviation Turbine Fuel – middle distillate
 - Kerosene-type (e.g. Jet A, A-1): C8-C18(C22)
 - Wide-cut (e.g. Jet B): C5-C16
- Strict upper limits on total aromatics (25%) and total di-aromatics (3%)
- No real "PIONA" analysis (matrix too complex)



FIA – Overview and Issues

- ASTM D1319-14 – Standard Test Method for Hydrocarbon Types in Liquid Petroleum Products by Fluorescent Indicator Adsorption
 - Silica gel fractionation not always clean cut between saturates, olefins, and aromatics
 - Performance issues with newer dye batches
 - Sole manufacturer of the dye no longer exists
 - FIA dye more valuable than gold???




Issues with ASTM D1840 Method Wavelength



- ASTM D1840 must use 285 nm for all absorbance measurements to eliminate the interference from mono-aromatics
- The relative absorbance for di-aromatics at <250 nm is greater than 10x the absorbance at 285 nm!
- The lower absolute response at 285 nm can lead to a large error range

4. Summary of Test Method

4.1 The total concentration of naphthalenes in jet fuels is determined by measurement of the absorbance at 285 nm of a solution of the fuel at known concentration.



1

Naphthalene

1-Methylnaphthalene

1,2-Dimethylnaphthalene

1,6,7-Trimethylnaphthalene

1,4,6,7-Tetramethylnaphthalene


Absorbance

0.8
0.6
0.4
0.2
0

125 145 165 185 205 225 245 265 285 305 325

Wavelength (nm)

0.08




Interferences in ASTM D1840 Di-Aromatic Determination

6. Interferences

6.1 Interferences add to the apparent naphthalene content. Phenanthrenes, dibenzothiophenes, biphenyls, benzothiophenes, and anthracenes interfere if present. The end point limitation of 315 °C will minimize this interference except for benzothiophenes and biphenyls. The contribution to measured naphthalene content by the presence of 1 % of such interfering compounds can be estimated from [Table 1](#).

TABLE 1 Interfering Compounds

Type of Interfering Compound	Error in Percentage of Naphthalenes Caused by 1 % Interfering Compound
Phenanthrenes	2
Dibenzothiophenes	2
Biphenyls	1
Benzothiophenes	0.6
Anthracenes	0.1



Caveats for Analyzing Jet Fuel using GC-VUV

- Olefins typically not present in jet fuel
 - However, their absorbance spectra have response in both the saturate (125-160 nm) and aromatic (170-200 nm) regions
 - Any reported values are likely S/A coelutions; therefore all olefin spectra removed from spectra library
- Classification of saturates nearly impossible at higher C#
 - Number of isomers increases dramatically with C#
 - At higher C# ($>C_{10}$), absorbance spectra of isoparaffins and naphthenes are too similar to distinguish in a complex sample

C#	P/I Isomers
9	35
12	355
15	4347
18	60523

Why Can We Analyze Jet Fuel So Quickly?

- The requisite spectral data is simplified
 - At higher C# ($>C_{10}$), absorbance spectra of isoparaffins and naphthenes are too similar to distinguish in a complex sample → ***all saturates quantified together***
 - Mono- and di-aromatics matched using representative absorbance spectra → ***specific compound identification not necessary***
- Distinct wavelength regions of maximum absorbance → ***easy spectral deconvolution***
 - Saturates (125-160 nm)
 - Mono-aromatics (175-205 nm)
 - Di-aromatics (205-240 nm)

GC-VUV Conditions for Jet Fuel Analysis

- Agilent 6890 GC
 - 1 μL injection
 - Inlet: 250°C, split 100:1
 - Column: Restek Rxi-1ms (30 m x 0.25 mm, 0.25- μm)
 - Flow: **2 mL/min helium**, constant flow
 - Oven: **50°C, hold 0.1 min; 15°C/min to 260°C** (run time – **14.1 min**)
- VUV Analytics VGA-100
 - Makeup gas: 0.23 psi N_2
 - Flow cell and transfer line: 275°C
 - Acquisition rate: **7 spectra/sec**
 - Acquisition range: 125-240 nm




GC-VUV Conditions for Jet Fuel Analysis

- Agilent 6890 GC
 - 1 μL injection
 - Inlet: 250°C, split 100:1
 - Column: Restek Rxi-1ms (30 m x 0.25 mm, 0.25- μm)
 - Flow: **2 mL/min helium**, constant flow
 - Oven: **50°C hold 0.5 min; 15°C/min to 250°C** (run time – **14.1 min**)
- VUV Analytics VVA-100
 - Makeup gas: 0.23 psi N_2
 - Flow cell and transfer line: 275°C
 - Acquisition rate: **7 spectra/sec**
 - Acquisition range: 125-240 nm

**SAME CONFIGURATION
AS ASTM D8071**




ASTM 1803 Jet Fuel



Fractionation of Jet Fuels

Conditioning



Red = methylene chloride
Blue = hexane

Sample Loading



Yellow = sample
Blue = hexane

Aliphatic Elution




Blue = hexane

Aromatic Elution



Red = methylene chloride

Concentration




Blue = Fraction 1
(aliphatics) in hexane


Red = Fraction 2
(aromatics) in methylene
chloride









Jet Fuel – Aliphatic Fraction






Jet Fuel – Aromatic Fraction












Jet Fuel Analysis – Total Saturates





Jet Fuel Analysis – Total Mono-Aromatics




Jet Fuel Analysis – Total Di-Aromatics






Analyte Table


RRF Display	Analyte Name	Analyte Category	Fit Values	Spectra Ave.	Response (Sum, 4 Scans)	Ret. Index	^
Total Aromatics		Mono-Aromatic	0.7298	0.3713	0.2710	1400.0	
Total Di-Aromatics		Di-Aromatic	0.5872	0.2555	0.1500	1418.0	
Total Saturates		Saturate	1.0999	0.1984	0.2177	1390.0	



Total Aromatics



Total Di-Aromatics



Inter-Run Precision of Jet Fuel

Total Saturates				Total Aromatics				Total Di-Aromatics			
Sample	Average	Std Dev	%RSD	Sample	Average	Std Dev	%RSD	Sample	Average	Std Dev	%RSD
AITF T277	82.2	0.03	0.03	AITF T277	17.8	0.03	0.15	AITF T277	1.88	0.003	0.14
ASTM 1703	80.9	0.02	0.02	ASTM 1703	19.1	0.02	0.08	ASTM 1703	0.87	0.004	0.48
ASTM 1711	87.6	0.01	0.02	ASTM 1711	12.4	0.01	0.12	ASTM 1711	2.11	0.006	0.26
CL18 3008	82.2	0.08	0.10	CL18 3008	17.8	0.08	0.45	CL18 3008	1.12	0.01	1.21
CL18 3011	84.4	0.02	0.03	CL18 3011	15.6	0.02	0.13	CL18 3011	0.10	0.003	3.34
CL18 3012	86.9	0.02	0.03	CL18 3012	13.1	0.02	0.17	CL18 3012	0.99	0.006	0.59

Conclusions



- PIONA analysis of gasoline
 - Fast and accurate, spectral deconvolution is very powerful
 - ASTM D8071 (IP# coming soon!), ILS complete
- Conjugated diolefins in gasoline
 - Can be analyzed as part of D8071; LODs as low as 0.01%
- VUV VHA
 - Prototype, but results match well with D8071
- Analysis of jet fuel
 - Data correlates well with existing methods, high precision
 - ILS complete, ASTM and IP #s coming very soon!