



V U V A N A L Y T I C S

S C I E N C E I N A N E W L I G H T

Verified Hydrocarbon Analysis™ (VHA™): Automating Complex Hydrocarbon Analysis on the VUV Analyzer Platform (GC-VUV)

Alex Hodgson, Applications Manager
VUV Analytics, Inc



Overview

- Who is VUV Analytics?
- Detailed Hydrocarbon Analysis (DHA)
- Verified Hydrocarbon Analysis™ (VHA)
 - Fast
 - Easy to use
 - More accurate results
- Case Study
- VUV Analyzer™ Platform for Fuels
- Summary






Who is VUV Analytics?



Since it's founding in 2009...

HEADQUARTERS IN CEDAR PARK, TX

Just outside of Austin, TX, VUV Analytics HQ is home to the global R&D, manufacturing and operations.



7 AWARDS & 20+ PATENTS

VUV Analytics has won over 7 Awards for the VGA Detector from multiple industries - process, oil & gas, R&D and cannabis. VUV holds over 20+ global patents on our technology and software. :

GLOBAL NETWORK

Even though VUV Analytics is based in the United States, their global partner network of distributors and 2 OEM partners gives them a worldwide presence, that is still growing.



OIL & GAS CUSTOMERS



The VUV Analyzer(TM) for Fuels continues to grow its customer base on a daily basis:
9 out of the top 13 Oil & Gas Companies
7 out of the top 10 Fuel Refineries
10+ Major Chemical Companies

PUBLICATIONS & METHODS

VUV Analytics has 48+ published scientific papers on GC-VUV technology. VUV Analytics has 2 approved ASTM methods and 1 working number.



LEARN MORE AT VUVANALYTICS.COM.



Detailed Hydrocarbon Analysis





Detailed Hydrocarbon Analysis

Traditional Approaches

| DHA Methods | | | | | | |
|---------------------|--|-------------|-------------|-------------|---|-------|
| Method | D6729 | D6730 | D6733 | Fast DHA | D7900 | D5134 |
| Sample Scope | Spark-ignition fuels Oxygenate blends | | | | Straight naphthas Reformate Alkylate Crude oil (No Gasolines and FCC Naphtha) | |
| Column (m) | 100m | 100m + 5m | 50m | 40m | 50m | |
| Conc. Range (m/m) | 0.01 - 30 % | 0.01 - 30 % | 0.01 - 15 % | 0.01 - 30 % | 0.01 - 30 % | |
| Max Olefin Conc. | 25% | 25% | 20% | 20% | 20% | |
| Max FBP | 225 °C | 225 °C | 225 °C | 225 °C | 270 °C | |
| Analysis Time (min) | 140 | 174 | 135 - 162 | 28 | 118 | |



Detailed Hydrocarbon Analysis

Traditional Approaches

| DHA Methods | | | | | | |
|---------------------|--|-------------|-------------|-------------|---|-------|
| Method | D6729 | D6730 | D6733 | Fast DHA | D7900 | D5134 |
| Sample Scope | Spark-ignition fuels Oxygenate blends | | | | Straight naphthas Reformate Alkylate Crude oil (No Gasolines and FCC Naphtha) | |
| Column (m) | 100m | 100m + 5m | 50m | 40m | 50m | |
| Conc. Range (m/m) | 0.01 - 30 % | 0.01 - 30 % | 0.01 - 15 % | 0.01 - 30 % | 0.01 - 30 % | |
| Max Olefin Conc. | 25% | 25% | 20% | 20% | 20% | |
| Max FBP | 225 °C | 225 °C | 225 °C | 225 °C | 270 °C | |
| Analysis Time (min) | 140 | 174 | 135 - 162 | 28 | 118 | |

No separation of
vital oxygenates
and toluene



Detailed Hydrocarbon Analysis

Traditional Approaches

| DHA Methods | | | | | | |
|---------------------|--|-------------|-------------|-------------|---|-------|
| Method | D6729 | D6730 | D6733 | Fast DHA | D7900 | D5134 |
| Sample Scope | Spark-ignition fuels Oxygenate blends | | | | Straight naphthas Reformate Alkylate Crude oil (No Gasolines and FCC Naphtha) | |
| Column (m) | 100m | 100m + 5m | 50m | 40m | 50m | |
| Conc. Range (m/m) | 0.01 - 30 % | 0.01 - 30 % | 0.01 - 15 % | 0.01 - 30 % | 0.01 - 30 % | |
| Max Olefin Conc. | 25% | 25% | 20% | 20% | 20% | |
| Max FBP | 225 °C | 225 °C | 225 °C | 225 °C | 270 °C | |
| Analysis Time (min) | 140 | 174 | 135 - 162 | 28 | 118 | |

No 1-methylnaphthalene
/ tridecane separation



Detailed Hydrocarbon Analysis

Traditional Approaches

| DHA Methods | | | | | | |
|---------------------|--|-------------|-------------|-------------|---|-------|
| Method | D6729 | D6730 | D6733 | Fast DHA | D7900 | D5134 |
| Sample Scope | Spark-ignition fuels Oxygenate blends | | | | Straight naphthas Reformate Alkylate Crude oil (No Gasolines and FCC Naphtha) | |
| Column (m) | 100m | 100m + 5m | 50m | 40m | 50m | |
| Conc. Range (m/m) | 0.01 - 30 % | 0.01 - 30 % | 0.01 - 15 % | 0.01 - 30 % | 0.01 - 30 % | |
| Max Olefin Conc. | 25% | 25% | 20% | 20% | 20% | |
| Max FBP | 225 °C | 225 °C | 225 °C | 225 °C | 270 °C | |
| Analysis Time (min) | 140 | 174 | 135 - 162 | 28 | 118 | |

Requires other methods
for benzene, toluene,
and oxygenate analysis



Detailed Hydrocarbon Analysis

Traditional Approaches

| DHA Methods | | | | | | |
|---------------------|--|-------------|-------------|-------------|---|-------|
| Method | D6729 | D6730 | D6733 | Fast DHA | D7900 | D5134 |
| Sample Scope | Spark-ignition fuels Oxygenate blends | | | | Straight naphthas Reformate Alkylate Crude oil (No Gasolines and FCC Naphtha) | |
| Column (m) | 100m | 100m + 5m | 50m | 40m | 50m | |
| Conc. Range (m/m) | 0.01 - 30 % | 0.01 - 30 % | 0.01 - 15 % | 0.01 - 30 % | 0.01 - 30 % | |
| Max Olefin Conc. | 25% | 25% | 20% | 20% | 20% | |
| Max FBP | 225 °C | 225 °C | 225 °C | 225 °C | 270 °C | |
| Analysis Time (min) | 140 | 174 | 135 - 162 | 28 | 118 | |

Only major peaks are identified after toluene



Detailed Hydrocarbon Analysis

Traditional Approaches

| DHA Methods | | | | | | |
|---------------------|--|-------------|-------------|-------------|---|-------|
| Method | D6729 | D6730 | D6733 | Fast DHA | D7900 | D5134 |
| Sample Scope | Spark-ignition fuels Oxygenate blends | | | | Straight naphthas Reformate Alkylate Crude oil (No Gasolines and FCC Naphtha) | |
| Column (m) | 100m | 100m + 5m | 50m | 40m | 50m | |
| Conc. Range (m/m) | 0.01 - 30 % | 0.01 - 30 % | 0.01 - 15 % | 0.01 - 30 % | 0.01 - 30 % | |
| Max Olefin Conc. | 25% | 25% | 20% | 20% | 20% | |
| Max FBP | 225 °C | 225 °C | 225 °C | 225 °C | 270 °C | |
| Analysis Time (min) | 140 | 174 | 135 - 162 | 28 | 118 | |

No separation of
oxygenates



Detailed Hydrocarbon Analysis

Traditional Approaches

| DHA Methods | | | | | | |
|---------------------|-------------|--|-------------|-------------|---|-------|
| Method | D6729 | D6730 | D6733 | Fast DHA | D7900 | D5134 |
| Sample Scope | | Spark-ignition fuels Oxygenate blends | | | Straight naphthas Reformate Alkylate Crude oil (No Gasolines and FCC Naphtha) | |
| Column (m) | 100m | 100m + 5m | 50m | 40m | 50m | |
| Conc. Range (m/m) | 0.01 - 30 % | 0.01 - 30 % | 0.01 - 15 % | 0.01 - 30 % | 0.01 - 30 % | |
| Max Olefin Conc. | 25% | 25% | 20% | 20% | 20% | |
| Max FBP | 225 °C | 225 °C | 225 °C | 225 °C | 270 °C | |
| Analysis Time (min) | 140 | 174 | 135 - 162 | 28 | 118 | |




Detailed Hydrocarbon Analysis

Traditional Approaches

| DHA Methods | |
|---------------------|--|
| Method | D6730 |
| Sample Scope | Spark-ignition fuels Oxygenate blends |
| Column (m) | 100m + 5m |
| Conc. Range (m/m) | 0.01 - 30 % |
| Max Olefin Conc. | 25% |
| Max FBP | 225 °C |
| Analysis Time (min) | 174 |

- **Uses Gas Chromatography (GC) with a Flame Ionization Detector (FID)**
 - GC provides high resolution separation of individual hydrocarbon components
 - FID (mass sensitive) detection is used to determine quantity of each component
- **Peak Identification is a combination of relative retention time databases and manual corrections**
- **Speciation of compounds in gasoline-range fuels, though not all compounds are identified by name.**
- **Potential areas of improvement**
 - Faster
 - Easier
 - More Accurate






Introduction to GC-VUV Spectroscopy





What is Vacuum Ultraviolet (VUV) Spectroscopy?

It's all about the light



- Works in a part of the electromagnetic spectrum that has previously been difficult to commercialize
- Characterized by very short wavelength (125 – 240 nm), high energy absorbance
- Nearly every compound absorbs in this region (except He, Ar, H₂)
- Compounds that absorb in this region have unique spectral fingerprints



What is Vacuum Ultraviolet (VUV) Spectroscopy?

The VGA Family of Detectors

- **Data confidence through spectral identification**

- Unique spectra = unambiguous compound identification
- Easily deconvolve co-eluting analytes
- Clear and easy isomer differentiation

- **Excellent sensitivity**

- Low picogram

- **Excellent temporal resolution**

- Up to 75Hz sampling

- **Predictable linear response**

- 1st principle detection reduces calibration burdens

- **Reliable & Easy to use**


- No vacuums pumps
- No need for baseline resolved peaks



| | VGA-100 | VGA-101 |
|-----------------|--------------|--------------|
| Max Temp | 300°C | 430°C |
| λ Range | 125 - 240 nm | 125 - 430 nm |
| Acq Rate | 75 Hz | 75 Hz |



How does a VGA detector work?




- GC analytes are passed to the VGA detector and excited by deuterium lamp
- Electronic transitions create unique spectral fingerprints
- Data is acquired in 3D (time, absorbance and wavelength) in VUVision™ Software for analysis
- Data is both qualitative and quantitative




Spectra Demonstrate Class-Based Similarities


Paraffins




Isoparaffins




Olefins



Naphthenes




Aromatics





Spectral Deconvolution of Coeluting Hydrocarbons


Acquired Spectrum
2-Methyl-1,3-butadiene
trans-2-Pentene
Pentane





Time Interval Deconvolution™ (TID)

- Quantitation method using VUV Analyze™ Software
- Chromatogram is divided into small time intervals (typically <0.05 min)
- For each time interval, compare measured spectrum against reference spectra in designated library, best analyte(s) fit determined
- Quickly determines total response per analyte
 - Converted to Mass % and Volume %





Verified Hydrocarbon Analysis™





Verified Hydrocarbon Analysis™

A Better Approach to DHA


Typical DHA versus VHA


| Method | ASTM D6730 (DHA) GC-FID | VHA VUV Analyzer™ Platform GC-VUV |
|------------------|----------------------------|---|
| Run Time (min.) | 174 | 50 |
| Column | Single 100m | Single 60m |
| Tuning Precolumn | Yes (5m) | no |
| Data Processing | Automated | Hands-off Analysis |
| Data Review | Manual | Automated |
| Verification | Retention Time | Spectral Validation |



Fast



Easy

Accurate

















0.46 AU







0.85 AU


















VHA™ – Deconvolution of Critical Pairs



- ASTM D6730
 - “Cut column until you have identified critical separations”
- NOT required with VHA thanks to critical pairs identification
- Fast
 - VHA enables compressed chromatography
- Easy
 - Deconvolution software automatically identifies compounds
- Accurate
 - Human guesswork is never a factor



VHA™ – Deconvolution of Critical Pairs





VUV Analyze: VHA

Step 1: Load the method

Step 2: Load a data file or an entire folder of files

File Status

Last Modification Date: 2019-09-06 10:29:34

Run File (*.db) or Directory

C:\Program Files\VUV Analytics\VUV Analyze\Demos\VHA\Sample219_101519

View Report


Navigate to File

Analyze

Step 3: Click "Analyze"

Step 4: Obtain Results

VHA™ vs DHA: Naphtha Streams



| Category | Naphtha A | | | Naphtha B | | |
|------------------|-------------|-------------|-------------|-------------|--------|-------------|
| | VHA 1 | VHA 2 | DHA | VHA 1 | VHA 2 | DHA |
| Paraffin | 21.0 | 20.8 | 19.1 | 23.2 | 23.4 | 22.0 |
| Isoparaffin | 31.4 | 30.8 | 27.6 | 31.5 | 31.3 | 25.4 |
| Olefin | 0.223 | 0.0775 | 0.861 | 0.210 | 0.0692 | 0.746 |
| Naphthene | 32.9 | 33.7 | 19.3 | 33.2 | | 23.0 |
| Aromatic | 14.5 | 14.6 | 15.2 | 14.5 | 14.6 | 13.3 |

| Category | Naphtha A | | | Naphtha B | | |
|-------------------|-------------|--------------|-------------|--------------|--------------|-------------|
| | VHA 1 | VHA 2 | DHA | VHA 1 | VHA 2 | DHA |
| Cyclopentane | 0.000 | 0.000 | 0.000 | 0.000 | 0.000 | 0.106 |
| MeCyclopentane | 0.000 | 0.000 | 0.000 | 0.000 | 1.46 | 1.38 |
| Ethylbenzene | 0.337 | 0.337 | 0.337 | 0.335 | 0.335 | 0.345 |
| o-Xylene | 1.31 | 1.31 | 1.31 | 1.31 | 1.32 | 1.45 |
| m-Xylene | ND | ND | ND | 0.047 | 0.074 | ND |
| p-Xylene | 3.16 | 3.16 | 3.11 | 4.60 | 4.67 | 4.75 |
| o-Dichlorobenzene | 1.48 | 1.50 | 1.41 | 1.90 | 1.92 | 1.79 |
| m-Dichlorobenzene | 1.51 | 1.52 | 1.89 | 0.700 | 0.702 | 0.948 |
| m-Xylene | 2.84 | 2.99 | 1.61 | 2.56 | 2.64 | 1.33 |
| p-Xylene | 1.01 | 0.935 | 1.58 | 0.835 | 0.794 | 1.27 |
| o-Xylene | 1.38 | 1.42 | 0.905 | 1.02 | 1.01 | 0.704 |

VHA Unidentified: <1.5%
DHA Unidentified: >15%



VHA™ vs DHA: European Oxygenates

Mass %

| Compound | Gasoline A | | | Gasoline B | | |
|--------------------------------|--------------|--------------|------------|--------------|--------------|---------------|
| | VHA 1 | VHA 2 | DHA | VHA 1 | VHA 2 | DHA |
| Cyclopentane | 0.359 | 0.334 | 0.359 | 0.348 | 0.350 | 0.358 |
| Methyl tert-butyl ether | 1.56 | 1.43 | ND* | 2.08 | 1.99 | 0.285* |
| Methylcyclopentane | 1.76 | 1.77 | 1.66 | 1.82 | 1.83 | 1.67 |
| Benzene | 0.812 | 0.804 | 0.759 | 0.785 | 0.775 | 0.750 |
| Cyclohexane | 1.03 | 1.02 | 1.11 | 1.02 | 0.967 | 1.09 |
| tert-Amyl methyl ether | 0.065 | 0.047 | ND | 0.089 | 0.071 | 0.062 |
| Isooctane | 0.033 | 0.058 | 0.065 | 0.032 | 0.072 | 0.077 |
| MeCyclohexane | 2.94 | 2.95 | 3.10 | 2.90 | 2.92 | 3.02 |
| Toluene | 10.6 | 10.6 | 10.3 | 10.1 | 9.96 | 9.89 |
| Ethylbenzene | 2.47 | 2.50 | 2.45 | 2.42 | 2.41 | 2.39 |
| m-Xylene | 5.58 | 5.99 | 5.42 | 5.53 | 5.83 | 5.30 |
| p-Xylene | 2.70 | 2.44 | 2.42 | 2.60 | 2.33 | 2.43 |
| o-Xylene | 3.13 | 3.19 | 3.05 | 3.09 | 3.09 | 3.01 |

***Misidentified as cis-4-Methyl-2-pentene in DHA method and data review**



VHA™ vs DHA: Aromatic Concentrates (aka PyGas)

| Mass % Category | PyGas A | | | PyGas B | | | PyGas C | | |
|--------------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| | VHA 1 | VHA 2 | DHA | VHA 1 | VHA 2 | DHA | VHA 1 | VHA 2 | DHA |
| Paraffin | 3.58 | 3.62 | 3.75 | 3.06 | 3.23 | 3.40 | 0.345 | 0.479 | 0.507 |
| Isoparaffin | 4.12 | 4.55 | 4.60 | 3.61 | 3.84 | 4.16 | 1.26 | 1.54 | 2.02 |
| Olefin | 15.0 | 15.6 | 16.0 | 23.3 | 24.1 | 21.0 | 23.6 | 24.8 | 18.0 |
| Naphthene | 6.69 | 7.25 | 5.56 | 4.06 | 4.42 | 5.05 | 0.469 | 0.507 | 0.487 |
| Aromatic | 70.6 | 69.0 | 66.2 | 66.0 | 64.4 | 60.9 | 74.4 | 72.6 | 60.3 |



VHA™ vs DHA: Aromatic Concentrates (aka PyGas)

| Mass % Compound | PyGas A | | | PyGas B | | | PyGas C | | |
|--------------------|---------|-------|-------|---------|-------|-------|---------|-------|-------|
| | VHA 1 | VHA 2 | DHA | VHA 1 | VHA 2 | DHA | VHA 1 | VHA 2 | DHA |
| Cyclopentane | 0.846 | 0.852 | 0.853 | 0.533 | 0.525 | 0.571 | 0.071 | 0.081 | 0.086 |
| MeCyclopentane | 2.40 | 2.37 | 2.24 | 1.84 | 1.79 | 1.70 | 0.065 | 0.077 | 0.086 |
| Benzene | 55.6 | 54.0 | 51.9 | 50.9 | 49.4 | 47.2 | 41.6 | 40.1 | 37.3 |
| Cyclohexane | 0.872 | 1.04 | 1.27 | 0.627 | 0.770 | 0.937 | ND | ND | 0.025 |
| Isooctane | 0.040 | 0.053 | 0.053 | ND | 0.012 | 0.025 | 0.044 | 0.017 | 0.009 |
| MeCyclohexane | 0.228 | 0.249 | 0.231 | 0.160 | 0.165 | 0.178 | ND | 0.011 | 0.020 |
| Toluene | 9.21 | 8.77 | 8.54 | 7.83 | 7.40 | 6.99 | 5.36 | 5.10 | 4.70 |
| Ethylbenzene | 1.06 | 0.999 | 0.976 | 0.666 | 0.669 | 0.608 | 0.958 | 0.936 | 0.838 |
| m-Xylene | 0.243 | 0.231 | 0.223 | 0.196 | 0.204 | 0.188 | 0.206 | 0.200 | 0.169 |
| p-Xylene | 0.132 | 0.118 | 0.108 | 0.096 | 0.104 | 0.095 | 0.098 | 0.087 | 0.083 |
| o-Xylene | 0.160 | 0.160 | 0.165 | 0.207 | 0.214 | 0.203 | 0.103 | 0.095 | 0.104 |



VHA™ vs DHA: Aromatic Concentrates (aka PyGas)

Mass %

| Compound | PyGas A | | | PyGas B | | | PyGas C | | |
|-----------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-----------|
| | VHA 1 | VHA 2 | DHA | VHA 1 | VHA 2 | DHA | VHA 1 | VHA 2 | DHA |
| 1-Pentene | 0.656 | 0.709 | 0.609 | 0.924 | 0.976 | 0.986 | 1.11 | 1.16 | 1.07 |
| n-Pentane | 0.610 | 0.655 | 0.644 | 0.258 | 0.270 | 0.313 | 0.097 | 0.095 | 0.102 |
| Isoprene | 1.03 | 1.08 | 0.812 | 0.703 | 0.722 | 0.612 | 0.371 | 0.392 | 0.311 |
| t-1,3-C5diene | 0.910 | 0.954 | 0.608 | 1.24 | 1.27 | 0.892 | 0.624 | 0.619 | 0.413 |
| Cyclopentadiene | 0.370 | 0.414 | 0.485 | 0.432 | 0.490 | 0.532 | 0.750 | 0.814 | 1.35 |
| c-1,3-C5diene | 0.526 | 0.548 | 0.459 | 0.702 | 0.714 | 0.609 | 0.304 | 0.310 | 0.231 |
| Cyclopentene | 0.767 | 0.785 | 0.866 | 0.763 | 0.760 | 0.900 | 4.04 | 4.04 | 4.53 |
| 2-Methylpentane | 1.05 | 1.04 | 0.946 | 0.961 | 0.939 | 0.929 | 0.174 | 0.190 | 0.181 |
| 3-Methylpentane | 1.04 | 1.02 | 1.04 | 0.963 | 0.941 | 1.01 | 0.105 | 0.103 | 0.112 |
| n-Hexane | 2.49 | 2.44 | 2.41 | 2.33 | 2.26 | 2.33 | 0.030 | 0.040 | 0.065 |
| EtCyPentane | 0.580 | 0.569 | 0.570 | 0.371 | 0.370 | 0.365 | ND | ND | 0.009 |
| Styrene | 1.93 | 1.83 | 1.61 | 3.09 | 2.92 | 2.47 | 3.86 | 3.57 | 3.184 |
| DCPD | 3.95 | 3.85 | 3.68 | 7.47 | 7.21 | 6.72 | 7.19 | 6.99 | ND |
| Indene | 0.445 | 0.548 | 0.585 | 0.651 | 0.780 | 0.744 | 1.04 | 1.20 | 1.09 |



Verified Hydrocarbon Analysis™

Key Benefits

- Fast
 - 50 minutes from start to finish
- Easy to use
 - Data processing and review is completely automated
- Accurate
 - Spectral verification and retention time

• Application Scope

- Finished Gasoline
- Oxygenates in gasoline
- Naphtha
- Reformate
- Alkylate
- PyGas
- FCC gasoline
- LPG

• Greatest Applicability

- Refinery streams
- VHA Templates
- On-line application



CASE STUDY

A major U.S. refinery suspected that their N2+A yield was low resulting in lower market prices for heavy naphtha

Closed Revenue
Leak (\$8 - \$16 million lost)




The VUV Analyzer™ Platform



The VUV Analyzer™ Platform

One Platform; Many Possibilities




ASTM D5599
GC-OFID

ASTM D5769
GC-MS

ASTM D3606
GC-TCD


ASTM D1319
FIA





The VUV Analyzer™ Platform

One Platform; Many Possibilities




Platform Advantages:

- One platform capable of running different methods
- No changes in hardware or setup
- Scalable from R&D to Production
- Fully automated applications
- Multi-method architecture
- Lower Total Cost of Analysis (CoA) compared to other techniques



Verified Hydrocarbon Analysis™ (VHA)



Summary

- Fast
 - VHA runtime is only 50 minutes
 - OVER 3X faster than traditional DHA
- Easy
 - Data processing and data review is fully automated with VUV Analyze™ Software
 - DHA requires skilled operator and manual data review
- Accurate
 - Data verification is performed using both retention time and spectral identification
 - DHA relies upon retention time alone



V U V A N A L Y T I C S

Alex Hodgson, VUV Analytics

Alex.Hodgson@vuvanalytics.com

512-333-0860

Follow me on LinkedIn: <https://www.linkedin.com/in/ahodgson91/>